

16S 454 Sequencing Protocol

HMP Consortium

Authors: Jumpstart Consortium Human Microbiome Project Data Generation Working Group

Version: 4.2.2

Effective Date: October 27, 2010

1 Abstract

2 Introduction

This SOP describes the 16S 454 Sequencing Protocol used by the HMP Consortium. The protocol describes the procedure for the clinical sample pilot study using barcoded primers for the 16S variable regions V1-3 and V3-5.

3 Requirements

3.1 Production Specifications

1. No more than one sample per 16S region per barcode will be sequenced per machine run.
2. All amplification pools will include one negative control (water) reaction. Centers should attempt to rotate the primer pair used for the negative control so to not always use the same primer pair.
3. All amplification pools will include one positive control reaction. Centers should attempt to rotate the primer pair used for the positive control so to not always use the same primer pair. The choice of positive control template is at the discretion of the sequencing center.
4. To reduce the potential for primer contamination, working stocks for all primer pairs should be tested for contamination by using each primer pair in an amplification reaction that contains no template. No amplified product should be observed when reactions are examined by electrophoresis on an agarose gel. If working stocks are stamped in batches across numerous aliquot plates, each plate should provide less than one weeks' worth of reactions (e.g. one to five uses). A sample plate from each stamped batch should be tested for contamination.
5. Sequencing centers will quantify the samples received using a fluorescent based assay and record this information.
6. Two attempts will be made to amplify each sample. The attempts will differ in the amount of template used:
 - a. Attempt 1: 2 μ l of template

16S 454 Sequencing Protocol

HMP Consortium

Authors: Jumpstart Consortium Human Microbiome Project Data Generation Working Group

Version: 4.2.2

Effective Date: October 27, 2010

- b. Attempt 2: 5 µl of template
 - c. After two attempts the sample can be considered a failure and no further work needs to be completed. However, it is up to the discretion of the center to make further attempts. It is also up to the discretion of the center as to whether or not to attempt sequencing of weakly amplifying samples
7. 5,000 reads should be attempted for each amplicon from a HMP donor sample
8. Amplicons that produce fewer than 3,000 reads passing QC (“good reads”) can be sequenced a second time to reach the deliverable of 3,000 reads passing QC. A center can choose to use the same amplicon, or produce a new amplicon for sequencing.
9. If the minimum number of reads (3,000 passing QC) from an amplicon has not been achieved after two sequencing attempts, no further sequencing needs to be completed.
10. Image capture and signal processing should use the most current version of the fragment processing software (V 2.3 as July 2010) as the default
11. Passing QC metrics for reads are:
 - a. > 300 nt (raw read)
 - b. Minimum of 300 Q20 bases
12. Technical replication between centers- 2% of HMP donor samples will be sequenced at two centers
13. A sample spreadsheet containing metadata and library construction information as set up by the DACC will be completed by all centers and will accompany the submission of the *.sff files.

3.2 Reagent Requirements

Material/Equipment	Vendor	Catalog Number
AccuPrime™ Taq DNA Polymerase High Fidelity	Invitrogen	12346-086
Forward and Reverse Primers premixed 96 well thermocycler plate clear adhesive plate seals	Operon	custom order
DNase/RNase free water	-	-
Thermo Cycler	-	-
Vortex	-	-
Pipettes	-	-

16S 454 Sequencing Protocol

HMP Consortium

Authors: Jumpstart Consortium Human Microbiome Project Data Generation Working Group

Version: 4.2.2

Effective Date: October 27, 2010

Aerosol resistant pipette tips MinElute PCR Purification Kit	Qiagen	28004
Ampure (SPRI) Beads (60mL kit) 1x low TE, pH 8.0	Agencourt	A29152
Quant-IT ds DNA Assay, high sensitivity	Invitrogen/Molecular Probes	Q33120

3.3 Additional Documentation

- Quant-iT ds DNA Assay protocol (manufacturer's specifications)
- SybrGreen Assay protocol (manufacturer's specifications)
- poolingCalculator.xls
- MinElute PCR Purification Kit Manual

4 Procedure

The PCR will be carried out using AccuPrime Taq High Fidelity. It is not necessary to setup this reaction on ice, however, it is recommended.

4.1 PCR Primer Setup

Set up of 10uM primer plates (combining barcoded A primer with non-barcoded B primer- see appendix below for primer & tag sequences):

- 1:10 dilution of the 100uM stocks:
- For each variable region, set up a working primer plate with 90ul of 1x low TE
- Add 5ul of each barcoded primer A from 100uM plate to corresponding well position in 10uM dilution plate.
- Add 5ul of the 100uM of corresponding B adapter to each well of the 96-well plate (final concentration 10uM primer pair).
 - o Mix by pipetting up and down.
- Working concentration of 4uM (2uM each primer)
 - o Dilute the 10uM primers 1:2.5 in 1x low TE (add 150ul of TE to each well of the 10uM primer stock plate and mix
 - o Primers can be stamped out into multiple single use primer plates and stored at -20°C until ready to use.

4.2 PCR Setup – Mastermix

16S 454 Sequencing Protocol

HMP Consortium

Authors: Jumpstart Consortium Human Microbiome Project Data Generation Working Group

Version: 4.2.2

Effective Date: October 27, 2010

4.2.1. MasterMix contains the following amounts per sample:

13.85uL	RNAase/DNAse free water
2uL	10X AccuPrime PCR Buffer II
0.15uL	AccuPrime Taq Hifi

16uL Total Volume of master mix

4.2.2. Multiply all the volumes above by the amount of reactions needed plus 10%

4.2.3. Combine reagents in a 2mL micro centrifuge tube and vortex to mix completely. If more than 100 reactions are needed a 15mL tube should be used.

4.2.4. Using an automated pipette transfer 16uL of master mix into individual wells in the 96 well reaction plate.

4.2.5. Cover plate and spin in a centrifuge at 2000rpm to collect sample at the bottom of the wells.

4.3 PCR Setup

4.3.1. Transfer 2uL diluted DNA sample into the respective reaction wells.

4.3.2. Transfer 2uL of barcoded primers from primer plate to corresponding wells in 96 well PCR plate.

4.3.3. Securely seal with clear adhesive plate seal and vortex plate vigorously.

4.3.4. Spin briefly at 2000 rpm in a centrifuge.

4.3.5. Place in thermo cycler and cycle as follows:

95°C	2 min	} 30 cycles
95°C	20 sec	
50 or 56°C*	30 sec	
72°C	5min	
4°C	forever	

*56°C for V3-1, 50°C for V5-3

4.3.6. Clean PCR products using Agencourt AmPure Beads

i. Use Agencourt protocol: 1.8x volume beads (36ul beads) – follow manufacturer’s specifications.

4.3.7. Elute beads with 25ul 1x low TE, pH 8.0 and transfer to new 96 well plate

4.4 PCR Gel Analysis

(E-gel alternative using 1ul of PCR product - faster) - we will know from the Quantification step below if we have product so this step is actually optional.

4.4.1. In a new reaction plate add 1uL PCR product to 1uL 6X loading dye.

4.4.2. Cover, vortex to mix, briefly centrifuge to collect sample at the bottom of the well.

16S 454 Sequencing Protocol

HMP Consortium

Authors: Jumpstart Consortium Human Microbiome Project Data Generation Working Group

Version: 4.2.2

Effective Date: October 27, 2010

- 4.4.3. Prepare a 1% agarose 1X TAE gel with EtBr.
- 4.4.4. Load samples and run approximately 1 hour at 100V.
- 4.4.5. Capture gel image on gel-doc and retain for analysis.

4.5 PCR Product Quantification

- 4.5.1 Quantify PCR product using SYBR-Green Quantification or Quant-IT ds DNA high sensitivity assay according to the manufacturer's specifications.

4.6 PCR Pooling

- 4.6.1 Using values from the SYBR Green or Quant-IT quantification, calculate pooling amounts using the poolingCalculator.xls or according to the following formula:

$$\text{Amount (uL) of each sample} = ((\text{vol}/2) * (\text{min})) / \text{sampleconc}$$

Where:

Vol = total volume of each sample

Min = concentration in ng/ul of the sample with the lowest concentration

Sampleconc = concentration in ng/uL of target sample

- 4.6.2 Pool samples using a minimum transfer volume of 1uL. If less than 1uL is called for, a dilution must be made. If using the poolingCalculator.xls this will be accounted for.
- 4.6.3 Using a Qiagen minElute column, purify the pool according to the manufacturer's protocol.

(The Broad normalizes by converting all concentrations to molecules/ul. Determine which sample has the lowest concentration and then dilute all other samples to the same concentration. Pool equal volume of each (5- 10ul) sample and then concentrate using a Qiagen MinElute column (elution with 30uL, 1x low TE, pH 8.0).

4.7 Sample Transfer for 454 Library Completion

- 4.7.1 Proceed directly to the qPCR library step.
- 4.7.2 *Optional:* Enter emPCR using 1/4 the recommended primer concentration to avoid too many molecules amplified on bead. *Higher primer concentrations may result in high signal intensities during run lead to higher mixed reads and shorter read lengths. Image software updates may reduce or eliminate this concern.*

5 Implementation

6 Discussion

Additional Information

16S 454 Sequencing Protocol HMP Consortium


Authors: Jumpstart Consortium Human Microbiome Project Data Generation Working Group

Version: 4.2.2

Effective Date: October 27, 2010

Broad institute primer sequences including tags

Purpose: In this approach, we will 454 barcode sequences designed internally by the Broad (Pablo Alvarez and Will Brockman) between the A adapter and primer specific sequence (see picture below). Barcoded primer sets have been tested by the Broad Institute.



454 Protocol Figure 1: Depiction of generalized workflow for 454 sequencing at HMP sequencing centers.

16S 454 Sequencing Protocol HMP Consortium

Authors: Jumpstart Consortium Human Microbiome Project Data Generation Working Group
Version: 4.2.2
Effective Date: October 27, 2010

Barcoded oligos for V3->V1 directional sequencing.		
Added the R specific primer sequence at 3' end of barcode on "A" adapter sequence		
Added the F specific primer sequence at the 3' end of the "B" adapter sequence		
		"B" adapter oligo sequence +27F (AGAGTTTGATCCTGGCTCAG)
		CCTATCCCTGTGTGCTTGGCAGTCTCAGAGAGTTTGATCCTGGCTCAG
Oligo name	Barcode	"A" adapter oligo sequence + barcode + 534R (ATTACCGGGCTGCTGG)
XLR_534R_v2bBar8L	CACGC	CCATCTCATCCCTGCGTGTCTCCGACTCAGCAGCATTACCGGGCTGCTGG
XLR_534R_v2bBar23L	CGCAAC	CCATCTCATCCCTGCGTGTCTCCGACTCAGCGCAACATTACCGGGCTGCTGG
XLR_534R_v2bBar174L	TGAAGC	CCATCTCATCCCTGCGTGTCTCCGACTCAGTGAAGCATTACCGGGCTGCTGG
XLR_534R_v2bBar602L	ACTTGC	CCATCTCATCCCTGCGTGTCTCCGACTCAGACTTGCATTACCGGGCTGCTGG
XLR_534R_v2bBar212L	TCACAC	CCATCTCATCCCTGCGTGTCTCCGACTCAGTTCACATTACCGGGCTGCTGG
XLR_534R_v2bBar25L	CGTGAC	CCATCTCATCCCTGCGTGTCTCCGACTCAGCGTGACATTACCGGGCTGCTGG
XLR_534R_v2bBar622L	ACGCGC	CCATCTCATCCCTGCGTGTCTCCGACTCAGACGCGCATTACCGGGCTGCTGG
XLR_534R_v2bBar72L	CCTCTC	CCATCTCATCCCTGCGTGTCTCCGACTCAGCCTCTCATTACCGGGCTGCTGG
XLR_534R_v2bBar600L	ACTCAC	CCATCTCATCCCTGCGTGTCTCCGACTCAGACTCACATTACCGGGCTGCTGG
XLR_534R_v2bBar559L	AGACAC	CCATCTCATCCCTGCGTGTCTCCGACTCAGAGACACATTACCGGGCTGCTGG
XLR_534R_v2bBar31L	CGACTC	CCATCTCATCCCTGCGTGTCTCCGACTCAGCGACTCATTACCGGGCTGCTGG
XLR_534R_v2bBar551L	AGCTTC	CCATCTCATCCCTGCGTGTCTCCGACTCAGAGCTTCATTACCGGGCTGCTGG
XLR_534R_v2bBar1149L	AAGCGCG	CCATCTCATCCCTGCGTGTCTCCGACTCAGAAGCGCGCATTACCGGGCTGCTGG
XLR_534R_v2bBar15L	CAAGAAC	CCATCTCATCCCTGCGTGTCTCCGACTCAGCAAGAACATTACCGGGCTGCTGG
XLR_534R_v2bBar556L	AGTTGGC	CCATCTCATCCCTGCGTGTCTCCGACTCAGAGTTGGCATTACCGGGCTGCTGG
XLR_534R_v2bBar144L	TATCAAC	CCATCTCATCCCTGCGTGTCTCCGACTCAGTATCAACATTACCGGGCTGCTGG
XLR_534R_v2bBar575L	AGGCGGC	CCATCTCATCCCTGCGTGTCTCCGACTCAGAGGCGGCATTACCGGGCTGCTGG
XLR_534R_v2bBar48L	CGGTATC	CCATCTCATCCCTGCGTGTCTCCGACTCAGCGGTATCATTACCGGGCTGCTGG
XLR_534R_v2bBar166L	TGACGAC	CCATCTCATCCCTGCGTGTCTCCGACTCAGTGACGACATTACCGGGCTGCTGG
XLR_534R_v2bBar613L	ACAAGGC	CCATCTCATCCCTGCGTGTCTCCGACTCAGACAAGGCATTACCGGGCTGCTGG
XLR_534R_v2bBar560L	AGACCTC	CCATCTCATCCCTGCGTGTCTCCGACTCAGAGACCTCATTACCGGGCTGCTGG
XLR_534R_v2bBar741L	ATACCAC	CCATCTCATCCCTGCGTGTCTCCGACTCAGATACCACATTACCGGGCTGCTGG
XLR_534R_v2bBar228L	TCGCGGC	CCATCTCATCCCTGCGTGTCTCCGACTCAGTTCGCGGCATTACCGGGCTGCTGG
XLR_534R_v2bBar807L	ATCTTAC	CCATCTCATCCCTGCGTGTCTCCGACTCAGATCTTACATTACCGGGCTGCTGG
XLR_534R_v2bBar1273L	AACCAGC	CCATCTCATCCCTGCGTGTCTCCGACTCAGAACCAGCATTACCGGGCTGCTGG
XLR_534R_v2bBar441L	TTCGAGC	CCATCTCATCCCTGCGTGTCTCCGACTCAGTTCGAGCATTACCGGGCTGCTGG
XLR_534R_v2bBar1174L	AAGGTGC	CCATCTCATCCCTGCGTGTCTCCGACTCAGAAGGTGCATTACCGGGCTGCTGG
XLR_534R_v2bBar209L	TCTTGGC	CCATCTCATCCCTGCGTGTCTCCGACTCAGTCTTGGCATTACCGGGCTGCTGG
XLR_534R_v2bBar153L	TAATCTC	CCATCTCATCCCTGCGTGTCTCCGACTCAGTAATCTCATTACCGGGCTGCTGG
XLR_534R_v2bBar213L	TCACCTC	CCATCTCATCCCTGCGTGTCTCCGACTCAGTTCACCTCATTACCGGGCTGCTGG
XLR_534R_v2bBar298L	TCCGCTC	CCATCTCATCCCTGCGTGTCTCCGACTCAGTCCGCTCATTACCGGGCTGCTGG
XLR_534R_v2bBar146L	TATTGAC	CCATCTCATCCCTGCGTGTCTCCGACTCAGTATTGACATTACCGGGCTGCTGG
XLR_534R_v2bBar554L	AGTCGAC	CCATCTCATCCCTGCGTGTCTCCGACTCAGAGTCGACATTACCGGGCTGCTGG
XLR_534R_v2bBar646L	ACGGCTC	CCATCTCATCCCTGCGTGTCTCCGACTCAGACGGCTCATTACCGGGCTGCTGG
XLR_534R_v2bBar158L	TGCGTTC	CCATCTCATCCCTGCGTGTCTCCGACTCAGTTCGCTTCATTACCGGGCTGCTGG
XLR_534R_v2bBar207L	TCTCGAC	CCATCTCATCCCTGCGTGTCTCCGACTCAGTCTCGACATTACCGGGCTGCTGG
XLR_534R_v2bBar77L	CCAGGAC	CCATCTCATCCCTGCGTGTCTCCGACTCAGCCAGGACATTACCGGGCTGCTGG
XLR_534R_v2bBar601L	ACTCCTC	CCATCTCATCCCTGCGTGTCTCCGACTCAGACTCCTCATTACCGGGCTGCTGG
XLR_534R_v2bBar481L	TTCTTGC	CCATCTCATCCCTGCGTGTCTCCGACTCAGTTCCTGCATTACCGGGCTGCTGG
XLR_534R_v2bBar419L	TTCATAC	CCATCTCATCCCTGCGTGTCTCCGACTCAGTTCATACATTACCGGGCTGCTGG
XLR_534R_v2bBar26L	CGTCGTC	CCATCTCATCCCTGCGTGTCTCCGACTCAGCGTCGTCATTACCGGGCTGCTGG
XLR_534R_v2bBar1172L	AAGGCAC	CCATCTCATCCCTGCGTGTCTCCGACTCAGAAGGCACATTACCGGGCTGCTGG
XLR_534R_v2bBar1210L	AACAACTC	CCATCTCATCCCTGCGTGTCTCCGACTCAGAACAACATTACCGGGCTGCTGG
XLR_534R_v2bBar606L	ACACGGAC	CCATCTCATCCCTGCGTGTCTCCGACTCAGACACGGACATTACCGGGCTGCTGG
XLR_534R_v2bBar159L	TGCCGAAC	CCATCTCATCCCTGCGTGTCTCCGACTCAGTTCGCCAATATTACCGGGCTGCTGG
XLR_534R_v2bBar147L	TATTCGTC	CCATCTCATCCCTGCGTGTCTCCGACTCAGTATTTCGTCATTACCGGGCTGCTGG
XLR_534R_v2bBar141L	TAGGAATC	CCATCTCATCCCTGCGTGTCTCCGACTCAGTAGGAATCATTACCGGGCTGCTGG
XLR_534R_v2bBar119L	CCGGCCAC	CCATCTCATCCCTGCGTGTCTCCGACTCAGCCGGCCACATTACCGGGCTGCTGG
XLR_534R_v2bBar1379L	AATGGTAC	CCATCTCATCCCTGCGTGTCTCCGACTCAGAATGGTACATTACCGGGCTGCTGG
XLR_534R_v2bBar208L	TCTCCGTC	CCATCTCATCCCTGCGTGTCTCCGACTCAGTCTCCGTCATTACCGGGCTGCTGG
XLR_534R_v2bBar1267L	AACCTGGC	CCATCTCATCCCTGCGTGTCTCCGACTCAGAACCTGGCATTACCGGGCTGCTGG
XLR_534R_v2bBar637L	ACGAAGTC	CCATCTCATCCCTGCGTGTCTCCGACTCAGACGAAGTCATTACCGGGCTGCTGG
XLR_534R_v2bBar435L	TTCTGGC	CCATCTCATCCCTGCGTGTCTCCGACTCAGTTCCTGGCATTACCGGGCTGCTGG
XLR_534R_v2bBar1202L	AACACAAC	CCATCTCATCCCTGCGTGTCTCCGACTCAGAACAACAATTACCGGGCTGCTGG
XLR_534R_v2bBar413L	TTCTTGC	CCATCTCATCCCTGCGTGTCTCCGACTCAGTTCCTTGCATTACCGGGCTGCTGG
XLR_534R_v2bBar289L	TCCAAGTC	CCATCTCATCCCTGCGTGTCTCCGACTCAGTTCGAAGTCATTACCGGGCTGCTGG
XLR_534R_v2bBar433L	TTCGCGAC	CCATCTCATCCCTGCGTGTCTCCGACTCAGTTCGCGACATTACCGGGCTGCTGG
XLR_534R_v2bBar121L	CCGGTCGC	CCATCTCATCCCTGCGTGTCTCCGACTCAGCCGGTCGCATTACCGGGCTGCTGG

16S 454 Sequencing Protocol HMP Consortium

Authors: Jumpstart Consortium Human Microbiome Project Data Generation Working Group

Version: 4.2.2

Effective Date: October 27, 2010

XLR_534R_v2bBar669L	ACCTGAAC	CCATCTCATCCCTGCGTGTCTCCGACTCAGACCTGAACATTACCGCGGCTGCTGG
XLR_534R_v2bBar1156L	AAGAGTTC	CCATCTCATCCCTGCGTGTCTCCGACTCAGAAGAGTTTATTACCGCGGCTGCTGG
XLR_534R_v2bBar370L	TTGACAAC	CCATCTCATCCCTGCGTGTCTCCGACTCAGTTGACAACATTACCGCGGCTGCTGG
XLR_534R_v2bBar281L	TCCGAGAC	CCATCTCATCCCTGCGTGTCTCCGACTCAGTCCAGAACATTACCGCGGCTGCTGG
XLR_534R_v2bBar49L	CGGTCTTC	CCATCTCATCCCTGCGTGTCTCCGACTCAGCGGTCTTATTACCGCGGCTGCTGG
XLR_534R_v2bBar1173L	AAGGCCTC	CCATCTCATCCCTGCGTGTCTCCGACTCAGAAGGCCTTATTACCGCGGCTGCTGG
XLR_534R_v2bBar599L	ACTAATTC	CCATCTCATCCCTGCGTGTCTCCGACTCAGACTAATTCATTACCGCGGCTGCTGG
XLR_534R_v2bBar167L	TGACCGTC	CCATCTCATCCCTGCGTGTCTCCGACTCAGTGACCGTATTACCGCGGCTGCTGG
XLR_534R_v2bBar161L	TGTCGGAC	CCATCTCATCCCTGCGTGTCTCCGACTCAGTGTGGACATTACCGCGGCTGCTGG
XLR_534R_v2bBar580L	AGGTGTGC	CCATCTCATCCCTGCGTGTCTCCGACTCAGAGGTTGTATTACCGCGGCTGCTGG
XLR_534R_v2bBar629L	ACGAGAAC	CCATCTCATCCCTGCGTGTCTCCGACTCAGACGAGAACATTACCGCGGCTGCTGG
XLR_534R_v2bBar184L	TGGTGAAC	CCATCTCATCCCTGCGTGTCTCCGACTCAGTGGTGAACATTACCGCGGCTGCTGG
XLR_534R_v2bBar233L	TCGTTGTC	CCATCTCATCCCTGCGTGTCTCCGACTCAGTCGTTGTCATTACCGCGGCTGCTGG
XLR_534R_v2bBar364L	TTGTGTTT	CCATCTCATCCCTGCGTGTCTCCGACTCAGTTGTGTTTATTACCGCGGCTGCTGG
XLR_534R_v2bBar78L	CCACGGTC	CCATCTCATCCCTGCGTGTCTCCGACTCAGCCACGGTATTACCGCGGCTGCTGG
XLR_534R_v2bBar393L	TTGGAGGC	CCATCTCATCCCTGCGTGTCTCCGACTCAGTTGGAGGCATTACCGCGGCTGCTGG
XLR_534R_v2bBar350L	TTATCGGC	CCATCTCATCCCTGCGTGTCTCCGACTCAGTTATCGGCATTACCGCGGCTGCTGG
XLR_534R_v2bBar1164L	AAGAAGAC	CCATCTCATCCCTGCGTGTCTCCGACTCAGAAGAGACATTACCGCGGCTGCTGG
XLR_534R_v2bBar1196L	AACTGTTC	CCATCTCATCCCTGCGTGTCTCCGACTCAGACTGTTTATTACCGCGGCTGCTGG
XLR_534R_v2bBar411L	TTCTCAAC	CCATCTCATCCCTGCGTGTCTCCGACTCAGTTCTCAACATTACCGCGGCTGCTGG
XLR_534R_v2bBar6L	TTCTCTTC	CCATCTCATCCCTGCGTGTCTCCGACTCAGTTCCTTATTACCGCGGCTGCTGG
XLR_534R_v2bBar1031L	ATTTCGTAC	CCATCTCATCCCTGCGTGTCTCCGACTCAGATTTCGTACATTACCGCGGCTGCTGG
XLR_534R_v2bBar76L	CCTTCCGC	CCATCTCATCCCTGCGTGTCTCCGACTCAGCCTTCCGCATTACCGCGGCTGCTGG
XLR_534R_v2bBar555L	AGTCCGTC	CCATCTCATCCCTGCGTGTCTCCGACTCAGAGTCCGTCATTACCGCGGCTGCTGG
XLR_534R_v2bBar378L	TTGAACTC	CCATCTCATCCCTGCGTGTCTCCGACTCAGTTGAACTATTACCGCGGCTGCTGG
XLR_534R_v2bBar1225L	AACGAGGC	CCATCTCATCCCTGCGTGTCTCCGACTCAGAACGAGGCATTACCGCGGCTGCTGG
XLR_534R_v2bBar99L	CCGTTTAC	CCATCTCATCCCTGCGTGTCTCCGACTCAGCCGTTTACATTACCGCGGCTGCTGG
XLR_534R_v2bBar236L	TCGAGGAAC	CCATCTCATCCCTGCGTGTCTCCGACTCAGTCGAGGAACATTACCGCGGCTGCTGG
XLR_534R_v2bBar731L	ACCGGAAGC	CCATCTCATCCCTGCGTGTCTCCGACTCAGACCGGAAGCATTACCGCGGCTGCTGG
XLR_534R_v2bBar628L	ACGTTCCAC	CCATCTCATCCCTGCGTGTCTCCGACTCAGACGTTCCACATTACCGCGGCTGCTGG
XLR_534R_v2bBar1250L	AACGGAGTC	CCATCTCATCCCTGCGTGTCTCCGACTCAGAACGGAGTCATTACCGCGGCTGCTGG
XLR_534R_v2bBar438L	TTCGTTATC	CCATCTCATCCCTGCGTGTCTCCGACTCAGTTTCGTTATATTACCGCGGCTGCTGG
XLR_534R_v2bBar693L	ACCGTAATC	CCATCTCATCCCTGCGTGTCTCCGACTCAGACCGTAATCATTACCGCGGCTGCTGG
XLR_534R_v2bBar672L	ACCTTGGTC	CCATCTCATCCCTGCGTGTCTCCGACTCAGACCTTGGTCATTACCGCGGCTGCTGG
XLR_534R_v2bBar355L	TTAAGATTTC	CCATCTCATCCCTGCGTGTCTCCGACTCAGTTAAGATTTCATTACCGCGGCTGCTGG
XLR_534R_v2bBar187L	TGGTTGGTC	CCATCTCATCCCTGCGTGTCTCCGACTCAGTGGTTGGTCATTACCGCGGCTGCTGG
XLR_534R_v2bBar162L	TGTCGGTC	CCATCTCATCCCTGCGTGTCTCCGACTCAGTGGTTCGGTCATTACCGCGGCTGCTGG
XLR_534R_v2bBar1292L	AACCGTGTTC	CCATCTCATCCCTGCGTGTCTCCGACTCAGAACCGTGTTCATTACCGCGGCTGCTGG
27F/534R_000	CGTGTGACTG	CCATCTCATCCCTGCGTGTCTCCGACTCAGCGTGTGACTGATTACCGCGGCTGCTGG
27F/534R_001	CAGATACGAC	CCATCTCATCCCTGCGTGTCTCCGACTCAGCAGATACGACATTACCGCGGCTGCTGG
27F/534R_002	AGCTCGAGCG	CCATCTCATCCCTGCGTGTCTCCGACTCAGAGCTCGAGCGATTACCGCGGCTGCTGG
27F/534R_003	CTATCGAGAG	CCATCTCATCCCTGCGTGTCTCCGACTCAGTATCGAGAGATTACCGCGGCTGCTGG
27F/534R_004	CTGACTATCG	CCATCTCATCCCTGCGTGTCTCCGACTCAGTACTATCGATTACCGCGGCTGCTGG
27F/534R_005	ATATATAGCG	CCATCTCATCCCTGCGTGTCTCCGACTCAGATATATAGCGATTACCGCGGCTGCTGG
27F/534R_006	CAGTACGATG	CCATCTCATCCCTGCGTGTCTCCGACTCAGCAGTACGATGATTACCGCGGCTGCTGG
27F/534R_007	ACTCGCTAGC	CCATCTCATCCCTGCGTGTCTCCGACTCAGACTCGCTAGCATTACCGCGGCTGCTGG
27F/534R_008	TGAGTCTATC	CCATCTCATCCCTGCGTGTCTCCGACTCAGTGTACTATATTACCGCGGCTGCTGG
27F/534R_009	TAGCACTACT	CCATCTCATCCCTGCGTGTCTCCGACTCAGTAGCACTACTATTACCGCGGCTGCTGG
27F/534R_010	AGCGTACGTG	CCATCTCATCCCTGCGTGTCTCCGACTCAGAGCGTACGTGATTACCGCGGCTGCTGG
27F/534R_011	ACTCGTGTAC	CCATCTCATCCCTGCGTGTCTCCGACTCAGACTCGTGTACATTACCGCGGCTGCTGG
27F/534R_012	TCTACAGTAG	CCATCTCATCCCTGCGTGTCTCCGACTCAGTCTACAGTAGATTACCGCGGCTGCTGG
27F/534R_013	ACTATACATC	CCATCTCATCCCTGCGTGTCTCCGACTCAGACTATACATCATTACCGCGGCTGCTGG
27F/534R_014	TGCGCGAGTG	CCATCTCATCCCTGCGTGTCTCCGACTCAGTGCAGGAGTATTACCGCGGCTGCTGG
27F/534R_015	TCGCACACGT	CCATCTCATCCCTGCGTGTCTCCGACTCAGTCGCACACGATTACCGCGGCTGCTGG
27F/534R_016	AGCTATATCG	CCATCTCATCCCTGCGTGTCTCCGACTCAGAGCTATATCGATTACCGCGGCTGCTGG
27F/534R_017	ACGATCGTAT	CCATCTCATCCCTGCGTGTCTCCGACTCAGACGATCGTATATTACCGCGGCTGCTGG
27F/534R_018	TGCATATACG	CCATCTCATCCCTGCGTGTCTCCGACTCAGTGCATATACGATTACCGCGGCTGCTGG
27F/534R_019	CGAGACACTG	CCATCTCATCCCTGCGTGTCTCCGACTCAGCGAGACACTGATTACCGCGGCTGCTGG
27F/534R_020	TGTGCGCTAG	CCATCTCATCCCTGCGTGTCTCCGACTCAGTGTGCGCTAGATTACCGCGGCTGCTGG
27F/534R_021	TCGTACACGG	CCATCTCATCCCTGCGTGTCTCCGACTCAGTGTACACGGATTACCGCGGCTGCTGG
27F/534R_022	CACCTACCTAG	CCATCTCATCCCTGCGTGTCTCCGACTCAGACTCAGTACTAGATTACCGCGGCTGCTGG
27F/534R_023	TGTACAGCTC	CCATCTCATCCCTGCGTGTCTCCGACTCAGTGTACAGCTCATTACCGCGGCTGCTGG
27F/534R_024	CTGTCTGACG	CCATCTCATCCCTGCGTGTCTCCGACTCAGTGTCTGACGATTACCGCGGCTGCTGG
27F/534R_025	CACACTCGCG	CCATCTCATCCCTGCGTGTCTCCGACTCAGCACACTCGCGATTACCGCGGCTGCTGG
27F/534R_026	CGCTCGTCTG	CCATCTCATCCCTGCGTGTCTCCGACTCAGCGCTCGTCTGATTACCGCGGCTGCTGG
27F/534R_027	AGCGACGTCT	CCATCTCATCCCTGCGTGTCTCCGACTCAGAGCGACGTCTATTACCGCGGCTGCTGG
27F/534R_028	CTCACGACGC	CCATCTCATCCCTGCGTGTCTCCGACTCAGCTCACGACGCATTACCGCGGCTGCTGG
27F/534R_029	ATGTCAGTCG	CCATCTCATCCCTGCGTGTCTCCGACTCAGATGTCAGTTCGATTACCGCGGCTGCTGG

16S 454 Sequencing Protocol HMP Consortium

Authors: Jumpstart Consortium Human Microbiome Project Data Generation Working Group

Version: 4.2.2

Effective Date: October 27, 2010

27F/534R_030	TCATAGACAC	CCATCTCATCCCTGCGTGTCTCCGACTCAGTCATAGACACATTACCGCGGCTGCTGG
27F/534R_031	ATGTACGTGT	CCATCTCATCCCTGCGTGTCTCCGACTCAGATGTACGTGATTACCGCGGCTGCTGG
27F/534R_032	ATAGCGTGAG	CCATCTCATCCCTGCGTGTCTCCGACTCAGATAGCGTGAGATTACCGCGGCTGCTGG
27F/534R_033	TCTGTAGTCT	CCATCTCATCCCTGCGTGTCTCCGACTCAGTCTGTAGTCTATTACCGCGGCTGCTGG
27F/534R_034	TGATATCGTC	CCATCTCATCCCTGCGTGTCTCCGACTCAGTGATATCGTCATTACCGCGGCTGCTGG
27F/534R_035	TCACTACATG	CCATCTCATCCCTGCGTGTCTCCGACTCAGTCACTACATGATTACCGCGGCTGCTGG
27F/534R_036	AGATACGCAG	CCATCTCATCCCTGCGTGTCTCCGACTCAGAGATACGCAGATTACCGCGGCTGCTGG
27F/534R_037	TATGACTGAG	CCATCTCATCCCTGCGTGTCTCCGACTCAGTATGACTGAGATTACCGCGGCTGCTGG
27F/534R_038	AGCTGACTAG	CCATCTCATCCCTGCGTGTCTCCGACTCAGAGCTGACTAGATTACCGCGGCTGCTGG
27F/534R_039	CGCTACGCGC	CCATCTCATCCCTGCGTGTCTCCGACTCAGCGCTACGCGCATTACCGCGGCTGCTGG
27F/534R_040	ACTGAGTGAG	CCATCTCATCCCTGCGTGTCTCCGACTCAGACTGAGTGAGATTACCGCGGCTGCTGG
27F/534R_041	AGACGCTACT	CCATCTCATCCCTGCGTGTCTCCGACTCAGAGACGCTACTATTACCGCGGCTGCTGG
27F/534R_042	TATCTAGACG	CCATCTCATCCCTGCGTGTCTCCGACTCAGTATCTAGACGATTACCGCGGCTGCTGG
27F/534R_043	TCGTAATATC	CCATCTCATCCCTGCGTGTCTCCGACTCAGTCTGTAATATTACCGCGGCTGCTGG
27F/534R_044	TACAGTGAGC	CCATCTCATCCCTGCGTGTCTCCGACTCAGTACAGTGAGCATTACCGCGGCTGCTGG
27F/534R_045	ATCGATAGAC	CCATCTCATCCCTGCGTGTCTCCGACTCAGATCGATAGACATTACCGCGGCTGCTGG
27F/534R_046	AGCAGAGACG	CCATCTCATCCCTGCGTGTCTCCGACTCAGAGCAGAGACGATTACCGCGGCTGCTGG
27F/534R_047	CGACGTGCGC	CCATCTCATCCCTGCGTGTCTCCGACTCAGCGACGTGCGCATTACCGCGGCTGCTGG
27F/534R_048	CACTCTATCG	CCATCTCATCCCTGCGTGTCTCCGACTCAGCACTCTATCGATTACCGCGGCTGCTGG
27F/534R_049	TGCTCAGACG	CCATCTCATCCCTGCGTGTCTCCGACTCAGTCTCAGACGATTACCGCGGCTGCTGG
27F/534R_050	ACGATGCTCG	CCATCTCATCCCTGCGTGTCTCCGACTCAGACGATGCTCGATTACCGCGGCTGCTGG
27F/534R_051	TCGTAGCAGC	CCATCTCATCCCTGCGTGTCTCCGACTCAGTCTAGCAGCATTACCGCGGCTGCTGG
27F/534R_052	TCGCGCATCG	CCATCTCATCCCTGCGTGTCTCCGACTCAGTCTCGCATCGATTACCGCGGCTGCTGG
27F/534R_053	TCGACGCTCT	CCATCTCATCCCTGCGTGTCTCCGACTCAGTCTGACGCTCTATTACCGCGGCTGCTGG
27F/534R_054	CGACGCACAG	CCATCTCATCCCTGCGTGTCTCCGACTCAGCGACGCACAGATTACCGCGGCTGCTGG
27F/534R_055	TGCGTAGACT	CCATCTCATCCCTGCGTGTCTCCGACTCAGTCTGCGTAGACTATTACCGCGGCTGCTGG
27F/534R_056	AGTGTACTGC	CCATCTCATCCCTGCGTGTCTCCGACTCAGAGTGTACTGCATTACCGCGGCTGCTGG
27F/534R_057	CTAGACTCAG	CCATCTCATCCCTGCGTGTCTCCGACTCAGCTAGACTCAGATTACCGCGGCTGCTGG
27F/534R_058	AGCGCTGTAG	CCATCTCATCCCTGCGTGTCTCCGACTCAGAGCGCTGTAGATTACCGCGGCTGCTGG
27F/534R_059	TCTCGAGAGT	CCATCTCATCCCTGCGTGTCTCCGACTCAGTCTCGAGAGATTACCGCGGCTGCTGG
27F/534R_060	CGAGTCCGAG	CCATCTCATCCCTGCGTGTCTCCGACTCAGCGAGTCCGAGATTACCGCGGCTGCTGG
27F/534R_061	TAGCTAGTAT	CCATCTCATCCCTGCGTGTCTCCGACTCAGTAGTATAGTATTACCGCGGCTGCTGG
27F/534R_062	AGAGTCCGCG	CCATCTCATCCCTGCGTGTCTCCGACTCAGAGAGTCCGCGATTACCGCGGCTGCTGG
27F/534R_063	CTCGTCACTG	CCATCTCATCCCTGCGTGTCTCCGACTCAGCTCGTCACTGATTACCGCGGCTGCTGG
27F/534R_064	AGTCTAGTCT	CCATCTCATCCCTGCGTGTCTCCGACTCAGAGTCTAGTCTATTACCGCGGCTGCTGG
27F/534R_065	TGTACTCACT	CCATCTCATCCCTGCGTGTCTCCGACTCAGTGTACTCACTATTACCGCGGCTGCTGG
27F/534R_066	CTATGTACAG	CCATCTCATCCCTGCGTGTCTCCGACTCAGCTATGTACAGATTACCGCGGCTGCTGG
27F/534R_067	TCGTGATAGT	CCATCTCATCCCTGCGTGTCTCCGACTCAGTCTGATAGTATTACCGCGGCTGCTGG
27F/534R_068	TGTGTACGAG	CCATCTCATCCCTGCGTGTCTCCGACTCAGTGTGTACGAGATTACCGCGGCTGCTGG
27F/534R_069	ATCTAGTCAC	CCATCTCATCCCTGCGTGTCTCCGACTCAGATCTAGTCACATTACCGCGGCTGCTGG
27F/534R_070	TATGAGAGTG	CCATCTCATCCCTGCGTGTCTCCGACTCAGTATGAGAGTGATTACCGCGGCTGCTGG
27F/534R_071	TACTGCTCAG	CCATCTCATCCCTGCGTGTCTCCGACTCAGTACTGCTCAGATTACCGCGGCTGCTGG
27F/534R_072	CTATACTACT	CCATCTCATCCCTGCGTGTCTCCGACTCAGCTATACTACTATTACCGCGGCTGCTGG
27F/534R_073	ACAGTGCTAC	CCATCTCATCCCTGCGTGTCTCCGACTCAGACAGTGCTACATTACCGCGGCTGCTGG
27F/534R_074	AGTATAGAGC	CCATCTCATCCCTGCGTGTCTCCGACTCAGAGTATAGAGCATTACCGCGGCTGCTGG
27F/534R_075	ACATCGCGGAG	CCATCTCATCCCTGCGTGTCTCCGACTCAGACATCGCGGAGATTACCGCGGCTGCTGG
27F/534R_076	ATGACGACTC	CCATCTCATCCCTGCGTGTCTCCGACTCAGATGACGACTCATTACCGCGGCTGCTGG
27F/534R_077	TGTATGTACT	CCATCTCATCCCTGCGTGTCTCCGACTCAGTGTATGTACTATTACCGCGGCTGCTGG
27F/534R_078	CGCGAGATAC	CCATCTCATCCCTGCGTGTCTCCGACTCAGCGGAGATACATTACCGCGGCTGCTGG
27F/534R_079	CTACAGTGTG	CCATCTCATCCCTGCGTGTCTCCGACTCAGCTACAGTGTGATTACCGCGGCTGCTGG
27F/534R_080	TATCACGATG	CCATCTCATCCCTGCGTGTCTCCGACTCAGTATCACGATGATTACCGCGGCTGCTGG
27F/534R_081	TGCTACGTCG	CCATCTCATCCCTGCGTGTCTCCGACTCAGTCTACGTCGATTACCGCGGCTGCTGG
27F/534R_082	CGTACGTGCG	CCATCTCATCCCTGCGTGTCTCCGACTCAGCGTACGTGCGATTACCGCGGCTGCTGG
27F/534R_083	TCAGCACTCG	CCATCTCATCCCTGCGTGTCTCCGACTCAGTCACTCGATTACCGCGGCTGCTGG
27F/534R_084	CGTGACTGCG	CCATCTCATCCCTGCGTGTCTCCGACTCAGCGTACTGCGATTACCGCGGCTGCTGG
27F/534R_085	ACTATAGTAC	CCATCTCATCCCTGCGTGTCTCCGACTCAGACTATAGTACATTACCGCGGCTGCTGG
27F/534R_086	ATATGTCTGT	CCATCTCATCCCTGCGTGTCTCCGACTCAGATATGTCTGTATTACCGCGGCTGCTGG
27F/534R_087	TACTAGATGT	CCATCTCATCCCTGCGTGTCTCCGACTCAGTACTAGATGTATTACCGCGGCTGCTGG
27F/534R_088	TGCTGTCTCT	CCATCTCATCCCTGCGTGTCTCCGACTCAGTGTCTCTCTATTACCGCGGCTGCTGG
27F/534R_089	CGTGACGATC	CCATCTCATCCCTGCGTGTCTCCGACTCAGCGTGACGATCATTACCGCGGCTGCTGG
27F/534R_090	CGCGTGTAC	CCATCTCATCCCTGCGTGTCTCCGACTCAGCGCGTGTACATTACCGCGGCTGCTGG
27F/534R_091	TCACGTATCT	CCATCTCATCCCTGCGTGTCTCCGACTCAGTACGTATCTATTACCGCGGCTGCTGG
27F/534R_092	TAGAGACTAG	CCATCTCATCCCTGCGTGTCTCCGACTCAGTAGAGACTAGATTACCGCGGCTGCTGG
27F/534R_093	TATGCGCGCG	CCATCTCATCCCTGCGTGTCTCCGACTCAGTATGCGCGCGATTACCGCGGCTGCTGG
27F/534R_094	CATATACACG	CCATCTCATCCCTGCGTGTCTCCGACTCAGCATATACACGATTACCGCGGCTGCTGG
27F/534R_095	TCGACTCGAT	CCATCTCATCCCTGCGTGTCTCCGACTCAGTCTCGATATTACCGCGGCTGCTGG
27F/534R_096	ACACAGTCGT	CCATCTCATCCCTGCGTGTCTCCGACTCAGACACAGTCGATTACCGCGGCTGCTGG
27F/534R_097	AGTACACGTC	CCATCTCATCCCTGCGTGTCTCCGACTCAGAGTACACGTCATTACCGCGGCTGCTGG

16S 454 Sequencing Protocol HMP Consortium

Authors: Jumpstart Consortium Human Microbiome Project Data Generation Working Group

Version: 4.2.2

Effective Date: October 27, 2010

27F/534R_098	TAGCGATGAC	CCATCTCATCCCTGCGTGTCTCCGACTCAGTAGCGATGACATTACCGCGGCTGCTGG
27F/534R_099	TGCGTATAGC	CCATCTCATCCCTGCGTGTCTCCGACTCAGTGCATATAGCATTACCGCGGCTGCTGG
27F/534R_100	CGACGCGATG	CCATCTCATCCCTGCGTGTCTCCGACTCAGCGACGCGATGATTACCGCGGCTGCTGG
27F/534R_101	ACGCACTGCG	CCATCTCATCCCTGCGTGTCTCCGACTCAGACGCACTGCGATTACCGCGGCTGCTGG
27F/534R_102	ACTGTGACTC	CCATCTCATCCCTGCGTGTCTCCGACTCAGACTGTGACTCATTACCGCGGCTGCTGG
27F/534R_103	TGATCGACAG	CCATCTCATCCCTGCGTGTCTCCGACTCAGTGTGACAGATTACCGCGGCTGCTGG
27F/534R_104	TAGTATCGAT	CCATCTCATCCCTGCGTGTCTCCGACTCAGTAGTATCGATATTACCGCGGCTGCTGG
27F/534R_105	TAGACGCATC	CCATCTCATCCCTGCGTGTCTCCGACTCAGTAGACGCATCATTACCGCGGCTGCTGG
27F/534R_106	TATCGATCTC	CCATCTCATCCCTGCGTGTCTCCGACTCAGTATCGATCTCATTACCGCGGCTGCTGG
27F/534R_107	TATCAGTCGT	CCATCTCATCCCTGCGTGTCTCCGACTCAGTATCAGTCGTATTACCGCGGCTGCTGG
27F/534R_108	ACAGCTATAG	CCATCTCATCCCTGCGTGTCTCCGACTCAGACAGCTATAGATTACCGCGGCTGCTGG
27F/534R_109	CACTCTCGAC	CCATCTCATCCCTGCGTGTCTCCGACTCAGCACTCTCGACATTACCGCGGCTGCTGG
27F/534R_110	AGCTACTCTG	CCATCTCATCCCTGCGTGTCTCCGACTCAGAGCTACTCTGATTACCGCGGCTGCTGG
27F/534R_111	ATACGAGAGC	CCATCTCATCCCTGCGTGTCTCCGACTCAGATACGAGAGCATTACCGCGGCTGCTGG
27F/534R_112	ACGTCGCAGT	CCATCTCATCCCTGCGTGTCTCCGACTCAGACGTCGCAGTATTACCGCGGCTGCTGG
27F/534R_113	ATGTGCTACT	CCATCTCATCCCTGCGTGTCTCCGACTCAGATGTGCTACTATTACCGCGGCTGCTGG
27F/534R_114	CATGTACGTC	CCATCTCATCCCTGCGTGTCTCCGACTCAGCATGTACGTCATTACCGCGGCTGCTGG
27F/534R_115	CACGCGTCTC	CCATCTCATCCCTGCGTGTCTCCGACTCAGCACGCGTCTCATTACCGCGGCTGCTGG
27F/534R_116	CGCTATCGAG	CCATCTCATCCCTGCGTGTCTCCGACTCAGCGCTATCGAGATTACCGCGGCTGCTGG
27F/534R_117	ACGACACGAG	CCATCTCATCCCTGCGTGTCTCCGACTCAGACGACACGAGATTACCGCGGCTGCTGG
27F/534R_118	TGCGCGTCGC	CCATCTCATCCCTGCGTGTCTCCGACTCAGTGCCTGCGATTACCGCGGCTGCTGG
27F/534R_119	TCAGCTCGTG	CCATCTCATCCCTGCGTGTCTCCGACTCAGTACGCTCGTGATTACCGCGGCTGCTGG
27F/534R_120	AGACGACTGT	CCATCTCATCCCTGCGTGTCTCCGACTCAGAGACGACTGTATTACCGCGGCTGCTGG
27F/534R_121	CACAGTATAC	CCATCTCATCCCTGCGTGTCTCCGACTCAGCACAGTATACATTACCGCGGCTGCTGG
27F/534R_122	ACGTCATCTG	CCATCTCATCCCTGCGTGTCTCCGACTCAGACGTCATCTGATTACCGCGGCTGCTGG
27F/534R_123	AGACTGTGAG	CCATCTCATCCCTGCGTGTCTCCGACTCAGAGACTGTGAGATTACCGCGGCTGCTGG
27F/534R_124	TACACATCAC	CCATCTCATCCCTGCGTGTCTCCGACTCAGTACACATCACATTACCGCGGCTGCTGG
27F/534R_125	ATAGCTCGAC	CCATCTCATCCCTGCGTGTCTCCGACTCAGATAGCTCGACATTACCGCGGCTGCTGG

16S 454 Sequencing Protocol HMP Consortium

Authors: Jumpstart Consortium Human Microbiome Project Data Generation Working Group
Version: 4.2.2
Effective Date: October 27, 2010

Barcoded oligos for V5->V3 directional sequencing.		
Added the R specific primer sequence at 3' end of barcode on "A" adapter sequence		
Added the F specific primer sequence at the 3' end of the "B" adapter sequence		
		"B" adapter oligo sequence +357F (CCTACGGGAGGCAGCAG)
		CCTATCCCTGTGTGCCTTGGCAGTCTCAGCCTACGGGAGGCAGCAG
Oligo name	Barcode	"A" adapter oligo sequence + barcode +926R (CCGTCGAATTCMTTTRAGT)
XLR_926R_v2bBar8L	CACGC	CCATCTCATCCCTGCGTGTCTCCGACTCAGCAGCCCGTCAATTCMTTTRAGT
XLR_926R_v2bBar23L	CGCAAC	CCATCTCATCCCTGCGTGTCTCCGACTCAGCGCAACCCGTCGAATTCMTTTRAGT
XLR_926R_v2bBar174L	TGAAGC	CCATCTCATCCCTGCGTGTCTCCGACTCAGTGAAGCCCGTCAATTCMTTTRAGT
XLR_926R_v2bBar602L	ACTTGC	CCATCTCATCCCTGCGTGTCTCCGACTCAGACTTGCCCGTCAATTCMTTTRAGT
XLR_926R_v2bBar212L	TCACAC	CCATCTCATCCCTGCGTGTCTCCGACTCAGTACACCCCGTCAATTCMTTTRAGT
XLR_926R_v2bBar25L	CGTGAC	CCATCTCATCCCTGCGTGTCTCCGACTCAGCGTGACCCGTCGAATTCMTTTRAGT
XLR_926R_v2bBar622L	ACGCGC	CCATCTCATCCCTGCGTGTCTCCGACTCAGACGCGCCCGTCAATTCMTTTRAGT
XLR_926R_v2bBar72L	CCTCTC	CCATCTCATCCCTGCGTGTCTCCGACTCAGCCTCTCCCGTCAATTCMTTTRAGT
XLR_926R_v2bBar600L	ACTCAC	CCATCTCATCCCTGCGTGTCTCCGACTCAGACTACCCCGTCAATTCMTTTRAGT
XLR_926R_v2bBar559L	AGACAC	CCATCTCATCCCTGCGTGTCTCCGACTCAGAGACCCCGTCAATTCMTTTRAGT
XLR_926R_v2bBar31L	CGACTC	CCATCTCATCCCTGCGTGTCTCCGACTCAGCGACTCCCGTCAATTCMTTTRAGT
XLR_926R_v2bBar551L	AGCTTC	CCATCTCATCCCTGCGTGTCTCCGACTCAGAGCTTCCCGTCAATTCMTTTRAGT
XLR_926R_v2bBar1149L	AAGCCGC	CCATCTCATCCCTGCGTGTCTCCGACTCAGAAGCCCGTCAATTCMTTTRAGT
XLR_926R_v2bBar15L	CAAGAAC	CCATCTCATCCCTGCGTGTCTCCGACTCAGCAAGAACCCTGTCGAATTCMTTTRAGT
XLR_926R_v2bBar556L	AGTTGGC	CCATCTCATCCCTGCGTGTCTCCGACTCAGAGTTGGCCCGTCAATTCMTTTRAGT
XLR_926R_v2bBar144L	TATCAAC	CCATCTCATCCCTGCGTGTCTCCGACTCAGTATCAACCCGTCGAATTCMTTTRAGT
XLR_926R_v2bBar575L	AGGCGGC	CCATCTCATCCCTGCGTGTCTCCGACTCAGAGGCGCCCGTCAATTCMTTTRAGT
XLR_926R_v2bBar48L	CGGTATC	CCATCTCATCCCTGCGTGTCTCCGACTCAGCGGTATCCCGTCAATTCMTTTRAGT
XLR_926R_v2bBar166L	TGACGAC	CCATCTCATCCCTGCGTGTCTCCGACTCAGTACGACCCCGTCAATTCMTTTRAGT
XLR_926R_v2bBar613L	ACAAGGC	CCATCTCATCCCTGCGTGTCTCCGACTCAGACAAGCCCGTCAATTCMTTTRAGT
XLR_926R_v2bBar560L	AGACCTC	CCATCTCATCCCTGCGTGTCTCCGACTCAGAGACCTCCCGTCAATTCMTTTRAGT
XLR_926R_v2bBar741L	ATACCAC	CCATCTCATCCCTGCGTGTCTCCGACTCAGATACCACCCGTCGAATTCMTTTRAGT
XLR_926R_v2bBar228L	TCGGGGC	CCATCTCATCCCTGCGTGTCTCCGACTCAGTTCGGGGCCCGTCAATTCMTTTRAGT
XLR_926R_v2bBar807L	ATCTTAC	CCATCTCATCCCTGCGTGTCTCCGACTCAGATCTTACCCTGTCGAATTCMTTTRAGT
XLR_926R_v2bBar1273L	AACCAGC	CCATCTCATCCCTGCGTGTCTCCGACTCAGAACCAGCCCGTCAATTCMTTTRAGT
XLR_926R_v2bBar441L	TTCGAGC	CCATCTCATCCCTGCGTGTCTCCGACTCAGTTCGAGCCCGTCAATTCMTTTRAGT
XLR_926R_v2bBar1174L	AAGGTGC	CCATCTCATCCCTGCGTGTCTCCGACTCAGAAGGTGCCCGTCAATTCMTTTRAGT
XLR_926R_v2bBar209L	TCTTGGC	CCATCTCATCCCTGCGTGTCTCCGACTCAGTCTTGGCCCGTCAATTCMTTTRAGT
XLR_926R_v2bBar153L	TAATCTC	CCATCTCATCCCTGCGTGTCTCCGACTCAGTAATCTCCCGTCAATTCMTTTRAGT
XLR_926R_v2bBar213L	TCACCTC	CCATCTCATCCCTGCGTGTCTCCGACTCAGTACCTCCCGTCAATTCMTTTRAGT
XLR_926R_v2bBar298L	TCCGCTC	CCATCTCATCCCTGCGTGTCTCCGACTCAGTCCGCTCCCGTCAATTCMTTTRAGT
XLR_926R_v2bBar146L	TATTGAC	CCATCTCATCCCTGCGTGTCTCCGACTCAGTATTGACCCCGTCAATTCMTTTRAGT
XLR_926R_v2bBar554L	AGTCGAC	CCATCTCATCCCTGCGTGTCTCCGACTCAGAGTCGACCCGTCGAATTCMTTTRAGT
XLR_926R_v2bBar646L	ACGGCTC	CCATCTCATCCCTGCGTGTCTCCGACTCAGACGGCTCCCGTCAATTCMTTTRAGT
XLR_926R_v2bBar158L	TGCGTTC	CCATCTCATCCCTGCGTGTCTCCGACTCAGTGCCTCCCGTCAATTCMTTTRAGT
XLR_926R_v2bBar207L	TCTCGAC	CCATCTCATCCCTGCGTGTCTCCGACTCAGTCTCGACCCCGTCAATTCMTTTRAGT
XLR_926R_v2bBar77L	CCAGGAC	CCATCTCATCCCTGCGTGTCTCCGACTCAGCCAGGACCCGTCGAATTCMTTTRAGT
XLR_926R_v2bBar601L	ACTCCTC	CCATCTCATCCCTGCGTGTCTCCGACTCAGACTCCTCCCGTCAATTCMTTTRAGT
XLR_926R_v2bBar481L	TTCTTGC	CCATCTCATCCCTGCGTGTCTCCGACTCAGTTCTTCCCGTCAATTCMTTTRAGT
XLR_926R_v2bBar419L	TTCATAC	CCATCTCATCCCTGCGTGTCTCCGACTCAGTTTCATACCCCGTCAATTCMTTTRAGT
XLR_926R_v2bBar26L	CGTCGTC	CCATCTCATCCCTGCGTGTCTCCGACTCAGCGTCGTCGTCGTCGAATTCMTTTRAGT
XLR_926R_v2bBar1172L	AAGGCAC	CCATCTCATCCCTGCGTGTCTCCGACTCAGAAGGCACCCGTCGAATTCMTTTRAGT
XLR_926R_v2bBar1210L	AACAATC	CCATCTCATCCCTGCGTGTCTCCGACTCAGAACAAATCCCGTCAATTCMTTTRAGT
XLR_926R_v2bBar606L	ACACGGAC	CCATCTCATCCCTGCGTGTCTCCGACTCAGACACGGACCCGTCGAATTCMTTTRAGT
XLR_926R_v2bBar159L	TGCCGAA	CCATCTCATCCCTGCGTGTCTCCGACTCAGTGCCGAAACCCGTCGAATTCMTTTRAGT
XLR_926R_v2bBar147L	TATTCGT	CCATCTCATCCCTGCGTGTCTCCGACTCAGTATTCGTCCCGTCAATTCMTTTRAGT
XLR_926R_v2bBar141L	TAGGAAT	CCATCTCATCCCTGCGTGTCTCCGACTCAGTAGGAATCCCGTCAATTCMTTTRAGT
XLR_926R_v2bBar119L	CCGGCCAC	CCATCTCATCCCTGCGTGTCTCCGACTCAGCCGGCCACCCGTCGAATTCMTTTRAGT
XLR_926R_v2bBar1379L	AATGGTAC	CCATCTCATCCCTGCGTGTCTCCGACTCAGAATGGTACCCGTCGAATTCMTTTRAGT
XLR_926R_v2bBar208L	TCTCCGTC	CCATCTCATCCCTGCGTGTCTCCGACTCAGTCTCCGTCGTCGTCGAATTCMTTTRAGT
XLR_926R_v2bBar1267L	AACCTGGC	CCATCTCATCCCTGCGTGTCTCCGACTCAGAACCCTGGCCCGTCAATTCMTTTRAGT
XLR_926R_v2bBar637L	ACGAAGTC	CCATCTCATCCCTGCGTGTCTCCGACTCAGACGAAGTCCCGTCAATTCMTTTRAGT
XLR_926R_v2bBar435L	TTCTGGC	CCATCTCATCCCTGCGTGTCTCCGACTCAGTTCTTCCCGTCAATTCMTTTRAGT
XLR_926R_v2bBar1202L	AACACAAC	CCATCTCATCCCTGCGTGTCTCCGACTCAGAACACAACCCGTCGAATTCMTTTRAGT
XLR_926R_v2bBar413L	TTCTTGC	CCATCTCATCCCTGCGTGTCTCCGACTCAGTTCTTCCCGTCAATTCMTTTRAGT
XLR_926R_v2bBar289L	TCCAAGTC	CCATCTCATCCCTGCGTGTCTCCGACTCAGTCCAAGTCCCGTCAATTCMTTTRAGT
XLR_926R_v2bBar433L	TTCCGAC	CCATCTCATCCCTGCGTGTCTCCGACTCAGTTCCCGACCCGTCGAATTCMTTTRAGT

16S 454 Sequencing Protocol HMP Consortium

Authors: Jumpstart Consortium Human Microbiome Project Data Generation Working Group
Version: 4.2.2
Effective Date: October 27, 2010

XLR 926R v2bBar121L	CCGGTCGC	CCATCTCATCCCTGCGTGTCTCCGACTCAGCCGGTCGCCCGTCAATTCMTTTRAGT
XLR 926R v2bBar669L	ACCTGAAC	CCATCTCATCCCTGCGTGTCTCCGACTCAGACCTGAACCCGTCAATTCMTTTRAGT
XLR 926R v2bBar1156L	AAGAGTTC	CCATCTCATCCCTGCGTGTCTCCGACTCAGAAGAGTTCGCCGTCAATTCMTTTRAGT
XLR 926R v2bBar370L	TTGACAAC	CCATCTCATCCCTGCGTGTCTCCGACTCAGTTGACAACCCGTCAATTCMTTTRAGT
XLR 926R v2bBar281L	TCCAGAAC	CCATCTCATCCCTGCGTGTCTCCGACTCAGTCCAGAACCCTCAATTCMTTTRAGT
XLR 926R v2bBar49L	CGGTCTTC	CCATCTCATCCCTGCGTGTCTCCGACTCAGCGGTCTCCCGTCAATTCMTTTRAGT
XLR 926R v2bBar1173L	AAGGCCTC	CCATCTCATCCCTGCGTGTCTCCGACTCAGAAGGCCTCCCGTCAATTCMTTTRAGT
XLR 926R v2bBar599L	ACTAATTC	CCATCTCATCCCTGCGTGTCTCCGACTCAGACTAATTCGCCGTCAATTCMTTTRAGT
XLR 926R v2bBar167L	TGACCGTC	CCATCTCATCCCTGCGTGTCTCCGACTCAGTGACCGTCCCGTCAATTCMTTTRAGT
XLR 926R v2bBar161L	TGTCGGAC	CCATCTCATCCCTGCGTGTCTCCGACTCAGTGTGACCGTCAATTCMTTTRAGT
XLR 926R v2bBar580L	AGGTTGTC	CCATCTCATCCCTGCGTGTCTCCGACTCAGAGTTGTCGCCGTCAATTCMTTTRAGT
XLR 926R v2bBar629L	ACGAGAAC	CCATCTCATCCCTGCGTGTCTCCGACTCAGACGAGAACCCTCAATTCMTTTRAGT
XLR 926R v2bBar184L	TGGTGAAC	CCATCTCATCCCTGCGTGTCTCCGACTCAGTGGTGAACCCGTCAATTCMTTTRAGT
XLR 926R v2bBar233L	TCGTGTGC	CCATCTCATCCCTGCGTGTCTCCGACTCAGTCGTTGCCGTCAATTCMTTTRAGT
XLR 926R v2bBar364L	TTGTGTTC	CCATCTCATCCCTGCGTGTCTCCGACTCAGTTGTGTGCCGTCAATTCMTTTRAGT
XLR 926R v2bBar78L	CCACGGTC	CCATCTCATCCCTGCGTGTCTCCGACTCAGCCACGGTCCCGTCAATTCMTTTRAGT
XLR 926R v2bBar393L	TTGGAGGC	CCATCTCATCCCTGCGTGTCTCCGACTCAGTTGGAGGCCGTCAATTCMTTTRAGT
XLR 926R v2bBar350L	TTATCGGC	CCATCTCATCCCTGCGTGTCTCCGACTCAGTTATCGGCCGTCAATTCMTTTRAGT
XLR 926R v2bBar1164L	AAGAAGAC	CCATCTCATCCCTGCGTGTCTCCGACTCAGAAGAAGACCCTCAATTCMTTTRAGT
XLR 926R v2bBar1196L	AACTGTTT	CCATCTCATCCCTGCGTGTCTCCGACTCAGAAGTTCGCCGTCAATTCMTTTRAGT
XLR 926R v2bBar411L	TTCTCAAC	CCATCTCATCCCTGCGTGTCTCCGACTCAGTTCTCAACCCTCAATTCMTTTRAGT
XLR 926R v2bBar6L	CTTCCTTC	CCATCTCATCCCTGCGTGTCTCCGACTCAGTTCCTTCGCCGTCAATTCMTTTRAGT
XLR 926R v2bBar1031L	ATTCTGTAC	CCATCTCATCCCTGCGTGTCTCCGACTCAGATTCGTACCCTCAATTCMTTTRAGT
XLR 926R v2bBar76L	CCTTCCGC	CCATCTCATCCCTGCGTGTCTCCGACTCAGCCTTCGCCGTCAATTCMTTTRAGT
XLR 926R v2bBar555L	AGTCCGTC	CCATCTCATCCCTGCGTGTCTCCGACTCAGAGTCCGTCCCGTCAATTCMTTTRAGT
XLR 926R v2bBar378L	TTGAAGTC	CCATCTCATCCCTGCGTGTCTCCGACTCAGTTGAAGTCCCGTCAATTCMTTTRAGT
XLR 926R v2bBar1225L	AACGAGGC	CCATCTCATCCCTGCGTGTCTCCGACTCAGAACGAGGCCCGTCAATTCMTTTRAGT
XLR 926R v2bBar99L	CCGTTTAC	CCATCTCATCCCTGCGTGTCTCCGACTCAGCCGTTACCCGTCAATTCMTTTRAGT
XLR 926R v2bBar236L	TCGAGGAC	CCATCTCATCCCTGCGTGTCTCCGACTCAGTCGAGGAAACCCTCAATTCMTTTRAGT
XLR 926R v2bBar731L	ACCGGAAGC	CCATCTCATCCCTGCGTGTCTCCGACTCAGACCGGAAGCCCGTCAATTCMTTTRAGT
XLR 926R v2bBar628L	ACGTTCCAC	CCATCTCATCCCTGCGTGTCTCCGACTCAGACGTTCCACCCTCAATTCMTTTRAGT
XLR 926R v2bBar1250L	AACGAGTC	CCATCTCATCCCTGCGTGTCTCCGACTCAGAACGAGTCCCGTCAATTCMTTTRAGT
XLR 926R v2bBar438L	TTCGTTATC	CCATCTCATCCCTGCGTGTCTCCGACTCAGTTTCGTTATCCCGTCAATTCMTTTRAGT
XLR 926R v2bBar693L	ACCGTAATC	CCATCTCATCCCTGCGTGTCTCCGACTCAGACCGTAATCCCGTCAATTCMTTTRAGT
XLR 926R v2bBar672L	ACCTTGGTC	CCATCTCATCCCTGCGTGTCTCCGACTCAGACCTTGGTCCCGTCAATTCMTTTRAGT
XLR 926R v2bBar355L	TTAAGATTC	CCATCTCATCCCTGCGTGTCTCCGACTCAGTTAAGATTCGCCGTCAATTCMTTTRAGT
XLR 926R v2bBar187L	TGGTTGGTC	CCATCTCATCCCTGCGTGTCTCCGACTCAGTGGTTGGTCCCGTCAATTCMTTTRAGT
XLR 926R v2bBar162L	TGTCGGTTC	CCATCTCATCCCTGCGTGTCTCCGACTCAGTGTCCGGTCCCGTCAATTCMTTTRAGT
XLR 926R v2bBar1292L	AACCGTGC	CCATCTCATCCCTGCGTGTCTCCGACTCAGAACCCTGTCGCCGTCAATTCMTTTRAGT
357F/926R 000	TCATAGACAG	CCATCTCATCCCTGCGTGTCTCCGACTCAGTCATAGACAGCCGTCAATTCMTTTRAGT
357F/926R 001	TATCACTACT	CCATCTCATCCCTGCGTGTCTCCGACTCAGTATCACTACGCCGTCAATTCMTTTRAGT
357F/926R 002	AGCGTCAGTA	CCATCTCATCCCTGCGTGTCTCCGACTCAGAGCGTCACTACCGTCAATTCMTTTRAGT
357F/926R 003	CTGTACGTAG	CCATCTCATCCCTGCGTGTCTCCGACTCAGTGTACGTAGCCGTCAATTCMTTTRAGT
357F/926R 004	AGTCTCTAGA	CCATCTCATCCCTGCGTGTCTCCGACTCAGAGTCTCTAGACCGTCAATTCMTTTRAGT
357F/926R 005	AGATACACAG	CCATCTCATCCCTGCGTGTCTCCGACTCAGAGATACACAGCCGTCAATTCMTTTRAGT
357F/926R 006	ACTCTAGTCT	CCATCTCATCCCTGCGTGTCTCCGACTCAGACTCTAGTCTCCGTCAATTCMTTTRAGT
357F/926R 007	AGTCAAGTGA	CCATCTCATCCCTGCGTGTCTCCGACTCAGAGTCAAGTGAACCGTCAATTCMTTTRAGT
357F/926R 008	CTACGTCTGT	CCATCTCATCCCTGCGTGTCTCCGACTCAGTACGTCTCTCCGTCAATTCMTTTRAGT
357F/926R 009	CGACTACGAG	CCATCTCATCCCTGCGTGTCTCCGACTCAGGACTACGAGCCGTCAATTCMTTTRAGT
357F/926R 010	TAGCACACTA	CCATCTCATCCCTGCGTGTCTCCGACTCAGTAGCACACTACCGTCAATTCMTTTRAGT
357F/926R 011	TACGAGTACA	CCATCTCATCCCTGCGTGTCTCCGACTCAGTACGAGTACACCGTCAATTCMTTTRAGT
357F/926R 012	TGCTACTGAG	CCATCTCATCCCTGCGTGTCTCCGACTCAGTGTACTGAGCCGTCAATTCMTTTRAGT
357F/926R 013	CACGATAGCG	CCATCTCATCCCTGCGTGTCTCCGACTCAGCAGATAGCGCCGTCAATTCMTTTRAGT
357F/926R 014	TATATCGACA	CCATCTCATCCCTGCGTGTCTCCGACTCAGTATATCGACCCGTCAATTCMTTTRAGT
357F/926R 015	TGTACTACAT	CCATCTCATCCCTGCGTGTCTCCGACTCAGTGTACTACATCCGTCAATTCMTTTRAGT
357F/926R 016	AGAGCGCGAG	CCATCTCATCCCTGCGTGTCTCCGACTCAGAGAGCGCGAGCCGTCAATTCMTTTRAGT
357F/926R 017	CGTAGATCGA	CCATCTCATCCCTGCGTGTCTCCGACTCAGCGTAGATCGACCGTCAATTCMTTTRAGT
357F/926R 018	TGATGACGCG	CCATCTCATCCCTGCGTGTCTCCGACTCAGTATGACGCGCCGTCAATTCMTTTRAGT
357F/926R 019	TCTCTCGAGA	CCATCTCATCCCTGCGTGTCTCCGACTCAGTCTCTCGAGCCGTCAATTCMTTTRAGT
357F/926R 020	TAGTGTAGCG	CCATCTCATCCCTGCGTGTCTCCGACTCAGTGTAGTGTAGCCCGTCAATTCMTTTRAGT
357F/926R 021	TCACGACGTA	CCATCTCATCCCTGCGTGTCTCCGACTCAGTACGACGACCGTCAATTCMTTTRAGT
357F/926R 022	TGTAGAGTAG	CCATCTCATCCCTGCGTGTCTCCGACTCAGTGTAGAGTAGCCGTCAATTCMTTTRAGT
357F/926R 023	TGCGTACTCA	CCATCTCATCCCTGCGTGTCTCCGACTCAGTGTACTACCGTCAATTCMTTTRAGT
357F/926R 024	ACGCACACGT	CCATCTCATCCCTGCGTGTCTCCGACTCAGACGACACGTCGCCGTCAATTCMTTTRAGT
357F/926R 025	TGAGTATGAG	CCATCTCATCCCTGCGTGTCTCCGACTCAGTGTATGAGCCGTCAATTCMTTTRAGT
357F/926R 026	TCTATACGCT	CCATCTCATCCCTGCGTGTCTCCGACTCAGTCTATACGTCGCCGTCAATTCMTTTRAGT
357F/926R 027	CAGTGAGACG	CCATCTCATCCCTGCGTGTCTCCGACTCAGCAGTGAGACGCCGTCAATTCMTTTRAGT

16S 454 Sequencing Protocol HMP Consortium

Authors: Jumpstart Consortium Human Microbiome Project Data Generation Working Group

Version: 4.2.2

Effective Date: October 27, 2010

357F/926R_028	CTAGTATGCG	CCATCTCATCCCTGCGTGTCTCCGACTCAGCTAGTATGCGCCGTC AATTTCMTTTRAGT
357F/926R_029	TCTACAGCGT	CCATCTCATCCCTGCGTGTCTCCGACTCAGTCTACAGCGTCCGTC AATTTCMTTTRAGT
357F/926R_030	ATCGCTAGTA	CCATCTCATCCCTGCGTGTCTCCGACTCAGATCGCTAGTACCCTCA AATTTCMTTTRAGT
357F/926R_031	AGCAGCTACG	CCATCTCATCCCTGCGTGTCTCCGACTCAGAGCAGCTACGCCGTC AATTTCMTTTRAGT
357F/926R_032	TCGCTATATA	CCATCTCATCCCTGCGTGTCTCCGACTCAGTCGCTATATACCCTCA AATTTCMTTTRAGT
357F/926R_033	TCGCTACGCG	CCATCTCATCCCTGCGTGTCTCCGACTCAGTCGCTACGCCGTC AATTTCMTTTRAGT
357F/926R_034	TGAGATCTCG	CCATCTCATCCCTGCGTGTCTCCGACTCAGTGAGATCTGCCGTC AATTTCMTTTRAGT
357F/926R_035	CGTGAGTGCG	CCATCTCATCCCTGCGTGTCTCCGACTCAGCGTGAGTGCGCCGTC AATTTCMTTTRAGT
357F/926R_036	TCGAGCACGT	CCATCTCATCCCTGCGTGTCTCCGACTCAGTCGAGCAGCTCCGTC AATTTCMTTTRAGT
357F/926R_037	AGACATATCG	CCATCTCATCCCTGCGTGTCTCCGACTCAGAGACATATGCCGTC AATTTCMTTTRAGT
357F/926R_038	TCTCGTGTGT	CCATCTCATCCCTGCGTGTCTCCGACTCAGTCTCGTGTGCCGTC AATTTCMTTTRAGT
357F/926R_039	ATATACGCGT	CCATCTCATCCCTGCGTGTCTCCGACTCAGATATACGCCGTCGTC AATTTCMTTTRAGT
357F/926R_040	ACGTCTCGCT	CCATCTCATCCCTGCGTGTCTCCGACTCAGACGCTCTCGCTCCGTC AATTTCMTTTRAGT
357F/926R_041	CGCGCTACGT	CCATCTCATCCCTGCGTGTCTCCGACTCAGCGCGCTACGTCGTC AATTTCMTTTRAGT
357F/926R_042	AGAGTCACGA	CCATCTCATCCCTGCGTGTCTCCGACTCAGAGAGTCAGACCGTC AATTTCMTTTRAGT
357F/926R_043	TACACGATGT	CCATCTCATCCCTGCGTGTCTCCGACTCAGTACAGATGTCCGTC AATTTCMTTTRAGT
357F/926R_044	CTACATCACG	CCATCTCATCCCTGCGTGTCTCCGACTCAGCTACATCAGCGGTC AATTTCMTTTRAGT
357F/926R_045	ATACTCTATG	CCATCTCATCCCTGCGTGTCTCCGACTCAGATACTCTATGCCGTC AATTTCMTTTRAGT
357F/926R_046	CTCACGCGAG	CCATCTCATCCCTGCGTGTCTCCGACTCAGCTCACGCCGAGCCGTC AATTTCMTTTRAGT
357F/926R_047	TCAGTCTCGA	CCATCTCATCCCTGCGTGTCTCCGACTCAGTCTCGTCTCGACCGTC AATTTCMTTTRAGT
357F/926R_048	AGAGTACGTG	CCATCTCATCCCTGCGTGTCTCCGACTCAGAGAGTACGTCGTC AATTTCMTTTRAGT
357F/926R_049	TGAGACGAGA	CCATCTCATCCCTGCGTGTCTCCGACTCAGTGAGACGAGACCGTC AATTTCMTTTRAGT
357F/926R_050	TATCTGTATA	CCATCTCATCCCTGCGTGTCTCCGACTCAGTATCTGTATACCCTCA AATTTCMTTTRAGT
357F/926R_051	TAGCGTGATG	CCATCTCATCCCTGCGTGTCTCCGACTCAGTAGCGTGATGCCGTC AATTTCMTTTRAGT
357F/926R_052	ACTGTATATG	CCATCTCATCCCTGCGTGTCTCCGACTCAGACTGTATATGCCGTC AATTTCMTTTRAGT
357F/926R_053	CATATAGACG	CCATCTCATCCCTGCGTGTCTCCGACTCAGCATATAGACGCCGTC AATTTCMTTTRAGT
357F/926R_054	ACGATAGACT	CCATCTCATCCCTGCGTGTCTCCGACTCAGACGATAGTCCGTC AATTTCMTTTRAGT
357F/926R_055	TGTCTGAGCA	CCATCTCATCCCTGCGTGTCTCCGACTCAGTGTCTGAGACCGTC AATTTCMTTTRAGT
357F/926R_056	AGTGACTAGT	CCATCTCATCCCTGCGTGTCTCCGACTCAGAGTGACTAGTCCGTC AATTTCMTTTRAGT
357F/926R_057	TCATCTAGT	CCATCTCATCCCTGCGTGTCTCCGACTCAGTCTATCTAGTCCGTC AATTTCMTTTRAGT
357F/926R_058	CTAGAGTGTG	CCATCTCATCCCTGCGTGTCTCCGACTCAGCTAGAGTGTGCCGTC AATTTCMTTTRAGT
357F/926R_059	TAGTACTGTG	CCATCTCATCCCTGCGTGTCTCCGACTCAGTAGTACTGTCCGTC AATTTCMTTTRAGT
357F/926R_060	ATACTGCGTA	CCATCTCATCCCTGCGTGTCTCCGACTCAGATACTGCGTACCCTCA AATTTCMTTTRAGT
357F/926R_061	TCAGCGTCTA	CCATCTCATCCCTGCGTGTCTCCGACTCAGTCAGCGTCTACCCTCA AATTTCMTTTRAGT
357F/926R_062	CGTACTCAGT	CCATCTCATCCCTGCGTGTCTCCGACTCAGCGTACTCAGTCCGTC AATTTCMTTTRAGT
357F/926R_063	TCGACGAGCT	CCATCTCATCCCTGCGTGTCTCCGACTCAGTCGACGAGCTCCGTC AATTTCMTTTRAGT
357F/926R_064	AGTCGACATG	CCATCTCATCCCTGCGTGTCTCCGACTCAGAGTCGACATGCCGTC AATTTCMTTTRAGT
357F/926R_065	CGTCAGCAGC	CCATCTCATCCCTGCGTGTCTCCGACTCAGCGTCAGCAGCCGTC AATTTCMTTTRAGT
357F/926R_066	TCGCTGATAG	CCATCTCATCCCTGCGTGTCTCCGACTCAGTCTGATAGCCGTC AATTTCMTTTRAGT
357F/926R_067	ACATACTGTG	CCATCTCATCCCTGCGTGTCTCCGACTCAGACATACTGTGCCGTC AATTTCMTTTRAGT
357F/926R_068	TCATCGAGTG	CCATCTCATCCCTGCGTGTCTCCGACTCAGTCTATCGAGTCCGTC AATTTCMTTTRAGT
357F/926R_069	ACTATATCTA	CCATCTCATCCCTGCGTGTCTCCGACTCAGACTATATCTACCCTCA AATTTCMTTTRAGT
357F/926R_070	ACTGTCTGTA	CCATCTCATCCCTGCGTGTCTCCGACTCAGACTGTCTGTACCCTCA AATTTCMTTTRAGT
357F/926R_071	CGTGTGCGAG	CCATCTCATCCCTGCGTGTCTCCGACTCAGCGTGTGCGAGCCGTC AATTTCMTTTRAGT
357F/926R_072	CTAGTGACAG	CCATCTCATCCCTGCGTGTCTCCGACTCAGCTAGTGCAGCCGTC AATTTCMTTTRAGT
357F/926R_073	TATACACGCG	CCATCTCATCCCTGCGTGTCTCCGACTCAGTATACAGCGCCGTC AATTTCMTTTRAGT
357F/926R_074	AGTCTAGTAG	CCATCTCATCCCTGCGTGTCTCCGACTCAGAGTCTAGTAGCCGTC AATTTCMTTTRAGT
357F/926R_075	CTCGTGACGT	CCATCTCATCCCTGCGTGTCTCCGACTCAGCTCGTGACGTCGTC AATTTCMTTTRAGT
357F/926R_076	TACACAGTCA	CCATCTCATCCCTGCGTGTCTCCGACTCAGTACACAGTACCCTCA AATTTCMTTTRAGT
357F/926R_077	CTACAGAGAG	CCATCTCATCCCTGCGTGTCTCCGACTCAGTACAGAGACCGTC AATTTCMTTTRAGT
357F/926R_078	TGACGTGACA	CCATCTCATCCCTGCGTGTCTCCGACTCAGTGACGTCGACCCGTC AATTTCMTTTRAGT
357F/926R_079	TCTGTGACAG	CCATCTCATCCCTGCGTGTCTCCGACTCAGTCTGTGACAGCCGTC AATTTCMTTTRAGT
357F/926R_080	TCTGCACTAG	CCATCTCATCCCTGCGTGTCTCCGACTCAGTCTGCACTAGCCGTC AATTTCMTTTRAGT
357F/926R_081	TGTACGACAG	CCATCTCATCCCTGCGTGTCTCCGACTCAGTGTACGACGCCGTC AATTTCMTTTRAGT
357F/926R_082	CTAGCTCGTA	CCATCTCATCCCTGCGTGTCTCCGACTCAGCTAGCTCGTACCCTCA AATTTCMTTTRAGT
357F/926R_083	ATACGCACGT	CCATCTCATCCCTGCGTGTCTCCGACTCAGATACGCAGCTCCGTC AATTTCMTTTRAGT
357F/926R_084	TGCTCTGCTAG	CCATCTCATCCCTGCGTGTCTCCGACTCAGTGTCTGTAGCCGTC AATTTCMTTTRAGT
357F/926R_085	AGTCAGGCGA	CCATCTCATCCCTGCGTGTCTCCGACTCAGAGTCGAGGACCGTC AATTTCMTTTRAGT
357F/926R_086	TATGACAGTG	CCATCTCATCCCTGCGTGTCTCCGACTCAGTATGACAGTCCGTC AATTTCMTTTRAGT
357F/926R_087	ACATAGTAGT	CCATCTCATCCCTGCGTGTCTCCGACTCAGACATAGTAGTCCGTC AATTTCMTTTRAGT
357F/926R_088	TATGATACTA	CCATCTCATCCCTGCGTGTCTCCGACTCAGTATGATACTACCCTCA AATTTCMTTTRAGT
357F/926R_089	TCGACGCATA	CCATCTCATCCCTGCGTGTCTCCGACTCAGTCGACGCATACCCTCA AATTTCMTTTRAGT
357F/926R_090	ACGCGAGATA	CCATCTCATCCCTGCGTGTCTCCGACTCAGACGCGAGATACCCTCA AATTTCMTTTRAGT
357F/926R_091	ACGATGATCG	CCATCTCATCCCTGCGTGTCTCCGACTCAGACGATGATGCCGTC AATTTCMTTTRAGT
357F/926R_092	ATCGTAGTGT	CCATCTCATCCCTGCGTGTCTCCGACTCAGATCTGATGTCGTC AATTTCMTTTRAGT
357F/926R_093	TATAGCGTCT	CCATCTCATCCCTGCGTGTCTCCGACTCAGTATAGCGTCTCCGTC AATTTCMTTTRAGT
357F/926R_094	ACTCTGTGAG	CCATCTCATCCCTGCGTGTCTCCGACTCAGACTCTGTGAGCCGTC AATTTCMTTTRAGT

16S 454 Sequencing Protocol HMP Consortium

Authors: Jumpstart Consortium Human Microbiome Project Data Generation Working Group

Version: 4.2.2

Effective Date: October 27, 2010

357F/926R_095	AGTAGCGTGT	CCATCTCATCCCTGCGTGTCTCCGACTCAGAGTAGCGTGTCCGTC AATTCMTTTRAGT
357F/926R_096	CGCGACGTGT	CCATCTCATCCCTGCGTGTCTCCGACTCAGCGCGACGTGTCCGTC AATTCMTTTRAGT
357F/926R_097	CTGTGTAGAG	CCATCTCATCCCTGCGTGTCTCCGACTCAGCTGTGTAGAGCCGTC AATTCMTTTRAGT
357F/926R_098	CTCTGTCTCG	CCATCTCATCCCTGCGTGTCTCCGACTCAGCTCTGTCTCGCCGTC AATTCMTTTRAGT
357F/926R_099	AGACGTCTCT	CCATCTCATCCCTGCGTGTCTCCGACTCAGAGACGTCTCTCCGTC AATTCMTTTRAGT
357F/926R_100	TCGAGAGTCG	CCATCTCATCCCTGCGTGTCTCCGACTCAGTCGAGAGTCGCCGTC AATTCMTTTRAGT
357F/926R_101	CTCTCGCGTA	CCATCTCATCCCTGCGTGTCTCCGACTCAGCTCTCGCGTACCCTGTC AATTCMTTTRAGT
357F/926R_102	ACGTGTAATA	CCATCTCATCCCTGCGTGTCTCCGACTCAGACGTGTAATAACCGTC AATTCMTTTRAGT
357F/926R_103	TGCTGCGTCT	CCATCTCATCCCTGCGTGTCTCCGACTCAGTGTGCTGCGTCCGTC AATTCMTTTRAGT
357F/926R_104	CATACTACTA	CCATCTCATCCCTGCGTGTCTCCGACTCAGCATACTACTACCGTC AATTCMTTTRAGT
357F/926R_105	AGTATCTCAG	CCATCTCATCCCTGCGTGTCTCCGACTCAGAGTATCTCAGCCGTC AATTCMTTTRAGT
357F/926R_106	TACTGCACAG	CCATCTCATCCCTGCGTGTCTCCGACTCAGTACTGCACAGCCGTC AATTCMTTTRAGT
357F/926R_107	CGCGCACGCG	CCATCTCATCCCTGCGTGTCTCCGACTCAGCGCGCACGCCGTC AATTCMTTTRAGT
357F/926R_108	TCACACTATA	CCATCTCATCCCTGCGTGTCTCCGACTCAGTCACACTATAACCGTC AATTCMTTTRAGT
357F/926R_109	TGACGCGCTA	CCATCTCATCCCTGCGTGTCTCCGACTCAGTGACGCGCTACCGTC AATTCMTTTRAGT
357F/926R_110	TGTACGTGTG	CCATCTCATCCCTGCGTGTCTCCGACTCAGTGTACGTGTGCCGTC AATTCMTTTRAGT
357F/926R_111	TCGTGATACA	CCATCTCATCCCTGCGTGTCTCCGACTCAGTCGTGATACACCGTC AATTCMTTTRAGT
357F/926R_112	ACACTATCAG	CCATCTCATCCCTGCGTGTCTCCGACTCAGACACTATCAGCCGTC AATTCMTTTRAGT
357F/926R_113	ATCGACGTCA	CCATCTCATCCCTGCGTGTCTCCGACTCAGATCGACGTCAACCGTC AATTCMTTTRAGT
357F/926R_114	TGATCTAGTA	CCATCTCATCCCTGCGTGTCTCCGACTCAGTGATCTAGTACCCTGTC AATTCMTTTRAGT
357F/926R_115	ATGACTGTCG	CCATCTCATCCCTGCGTGTCTCCGACTCAGATGACTGTGCCGTC AATTCMTTTRAGT
357F/926R_116	ATCGACAGAG	CCATCTCATCCCTGCGTGTCTCCGACTCAGATCGACAGAGCCGTC AATTCMTTTRAGT
357F/926R_117	AGTATATGTG	CCATCTCATCCCTGCGTGTCTCCGACTCAGAGTATATGTGCCGTC AATTCMTTTRAGT
357F/926R_118	AGACTGACAG	CCATCTCATCCCTGCGTGTCTCCGACTCAGAGACTGACAGCCGTC AATTCMTTTRAGT
357F/926R_119	AGCGCGTAGT	CCATCTCATCCCTGCGTGTCTCCGACTCAGAGCGCGTAGTCCGTC AATTCMTTTRAGT
357F/926R_120	AGACGATGTG	CCATCTCATCCCTGCGTGTCTCCGACTCAGAGACGATGTGCCGTC AATTCMTTTRAGT
357F/926R_121	TCAGTAGTGA	CCATCTCATCCCTGCGTGTCTCCGACTCAGTACAGTAGTACCGTC AATTCMTTTRAGT
357F/926R_122	TGCAGTCGTG	CCATCTCATCCCTGCGTGTCTCCGACTCAGTGCAGTCGTGCCGTC AATTCMTTTRAGT
357F/926R_123	TGCGTACATG	CCATCTCATCCCTGCGTGTCTCCGACTCAGTGCAGTACATGCCGTC AATTCMTTTRAGT
357F/926R_124	TCAGATGACG	CCATCTCATCCCTGCGTGTCTCCGACTCAGTGCAGATGACGCCGTC AATTCMTTTRAGT
357F/926R_125	AGTGTGAGTG	CCATCTCATCCCTGCGTGTCTCCGACTCAGAGTGTGAGTGCCGTC AATTCMTTTRAGT

7 Related Documents & References

8 Revision History

Version	Author/Reviewer	Date	Change Made
4.2		10/27/2010	Revised for Protocol
4.2.2		05/09/2012	Convert SOP to standard template